

# ADB-5Br-INACA



5-bromo-N-(1-carbamoyl-2,2-dimethyl-propyl)-1H-indazole-3-carboxamide

Formula: C<sub>14</sub>H<sub>17</sub>BrN<sub>4</sub>O<sub>2</sub> Formula weight: 353.21 Chemical Abstracts No.: *n. a.* Smiles code: CC(C)(C)C(NC(=O)c1n[nH]c2ccc(Br)cc12)C(=O)N InChi key: AJGASUCDTSLMNP-UHFFFAOYSA-N Other names:

#### Reference for the given name:

A. J. Potts, C. Cano, S. H. L. Thomas, S. L. Hill: Synthetic cannabinoid receptor agonists: classification and nomenclature, *Clinical Toxicology*, **58**(2), 89-92 (2020)

https://doi.org/10.1080/15563650.2019.1661425

The evidence was 4.58 grams off-white powder.

< No characteristic photo available >





## GC-MS

An Agilent 6890N Network GC system set up with Agilent HP-5MS (length: 30 m, diameter: 0.25 mm, film: 0.25 mm) coupled to an Agilent 5973 Network Mass Selective Detector (scan range m/z 35 – m/z 500) was used. The methanolic solution of the evidence was injected. Samples were subjected to electron ionization (EI) mode. GC-MS conditions: HP-5MS column was temperature programmed from 100  $^{\circ}$ C (which was held for 2 minutes) to 280  $^{\circ}$ C at 20  $^{\circ}$ C/min, 280  $^{\circ}$ C was held for 3 minutes, then to 315  $^{\circ}$ C at 25  $^{\circ}$ C/min, the temperature was stated at 315  $^{\circ}$ C for 12 minutes. The carrier gas was helium. Tribenzyl-amine was applied as an internal standard (locked to 10.8 minutes). Data handling was carried out with GC/MSD ChemStation software.



### GC-MS total ion current chromatogram

Agilent 6890N Network GC system set up with Agilent HP-5MS





## Mass spectrum at 15.60 min retention time

### Interpretation of the mass spectrum

Some peak pairs indicate bromine content in the molecule.

Agilent 6890N Network GC system set up with Agilent HP-5MS



## IR

The IR spectrum of the evidence was recorded on a Thermo SCIENTIFIC Nicolet iS5 FT-IR spectrometer equipped with an iD5 ATR accessory, in absorbance mode. The digital resolution is 4 cm<sup>-1</sup>. The spectrometer was controlled, and the data were processed using Omnic 9 software package. The spectrum was off-line imported into Bruker OPUS software, and the output below was performed by OPUS 7.5 software.



### IR spectrum of the evidence as received

Thermo SCIENTIFIC Nicolet iS5 FT-IR spectrometer



# NMR

The NMR spectra were recorded on a Bruker Avance Neo 400 NMR spectrometer operating at 9.4 Tesla magnetic field, equipped with Prodigy BBO-H&F-D-05 Z-gradient probe. The spectra were recorded at 25 °C in DMSO-*d*<sub>6</sub> solution. The spectrometer was controlled, and the data were processed using TopSpin 4.0 software package. Chemical shifts ( $\delta$ ) are given in parts per million unit, referenced to tetramethylsilane ( $\delta_{TMS}$  = 0.00 ppm). The determination of the structure was based on <sup>1</sup>H, zqs-easy-ROESY, as well as <sup>13</sup>C, multiplicity edited HSQC and HMBC spectra.

#### Interpretation of the NMR spectra



In DMSO-d<sub>6</sub> solution

Characteristic H-H steric proximities detected by zqs-easy-ROESY measurement



<sup>13</sup>C-NMR chemical shifts  $\delta$  [ppm]





Characteristic heteronuclear long-range couplings detected by HMBC measurement H ---- C \_\_\_







<sup>1</sup>H-NMR spectrum (overview and characteristic sections)





Bruker AVANCE NEO 400, CryoProbe Prodigy; solvent: DMSO-d<sub>6</sub>





Bruker AVANCE NEO 400, CryoProbe Prodigy; solvent: DMSO-d<sub>6</sub>





Bruker AVANCE NEO 400, CryoProbe Prodigy; solvent: DMSO-d<sub>6</sub>