


ANALYTICAL REPORT

Mephedrone (C₁₁H₁₅NO)

(RS)-2-Methylamino-1-(4-methylphenyl)propan-1-one

Remark – other active cpd. detected: **none**

Sample ID:	1000-12A
Sample description:	powder - white
Sample type:	RM-reference material
Comments:	Fluka Lot#BCBF2682V;
Date of entry:	2/2/2012

Substance identified-structure ¹ (base form)	
Systematic name:	(RS)-2-Methylamino-1-(4-methylphenyl)propan-1-one
Other names:	4-methylmethcathinone (4-MMC)
Formula (per base form)	C ₁₁ H ₁₅ NO
M _w (g/mol)	177,24
Salt form:	HCl
StdInChIKey	YELGFTGWJGBAQU-UHFFFAOYSA-N
Compound Class	Cathinones
Other active cpd. detected	none
Add.info (purity..)	99,90%

¹ Created by OPSIN free tool: <http://opsin.ch.cam.ac.uk/> DOI: 10.1021/ci100384d



Report updates

date	comments (explanation)

Supporting information

Analytical technique:	applied	remarks
GC-MS (EI ionization)	+	NFL GC-RT (min): 3,52 BP(1): 58; BP(2): 91,BP(3) :32,
FTIR-ATR	+	direct measurement
GC-IR (condensed phase)	/	

GC-MS (Agilent):

GC-method is RT locked to tetracosane (RT=9.53 min).

Injection volume 1 ml and split mode (1:50) .

Injector temperature: 280 °C.

Chromatographic separation

Column: HP1-MS (100% dimethylpolysiloxane), length 30 m, internal diameter 0.25 mm, film thickness 0.25 mm.

Carrier gas He: flow-rate 1.2 ml/min. GC oven program: 170 °C for 1 min, followed by heating up to 293 °C at a rate of 18 °C/min, hold for 6.1 min, then heating at 50 °C/min up to 325 °C and finally 2.8 min isothermal.

MSD source EI = 70 eV. GC-MS transfer line T= 235°C, source and quadropole temperatures 280°C and 180°C, respectively. Scan range m/z scan range: from 50 (40) to 550 amu.

FTIR-ATR (Perkin Elmer): scan range 4000-400 cm⁻¹; resolution 4cm⁻¹

GC- (MS)-IR condensed phase (GC-MS (Agilent) & IR (Spectra analyses-Danny) IR scan range 4000 to 700, resolution 4cm⁻¹

GC-method:

Injection volume 1 ml and split mode (1:5) .

Injector temperature: 280 °C.

Chromatographic separation

Column: HP1-MS (100% dimethylpolysiloxane), length 30 m, internal diameter 0.25 mm, film thickness 0.25 mm.

Carrier gas He: flow-rate 1.2 ml/min. GC oven program: 170 °C for 1 min, followed by heating up to 293 °C at a rate of 18 °C/min, hold for 6.1 min, then heating at 50 °C/min up to 325 °C and finally 2.8 min isothermal.

Split MS : IR : (1:9)


MSD source EI = 70 eV. GC-MS transfer line T= 235°C, source and quadropole temperatures 280°C and 180°C, respectively. Scan range m/z scan range: from 50 (40) to 550 amu.

IR (condensed phase): IR scan range 4000 to 700, resolution 4cm⁻¹


FIGURES OF SPECTRA

MS (EI)

Abundance



FTIR-ATR



Name	Description
Mephedrone HCl (4-MMC)_Fluka_Lot#BCBF2682V_ID 1000-12A	4-Methylmethcathinone HCl_ATR