


ANALYTICAL REPORT

Phenazepam-4-oxide (C₁₅H₁₀BrClN₂O₂)

7-bromo-5-(2-chlorophenyl)-4-hydroxy-3,4-dihydro-2H-1,4-benzodiazepin-2-one Phenazepam-4-oksido (Chiron name), Phenazepam-4-hydroxy

Remark – other active cpd. detected: **phenazepam**

Sample ID:	1236-15
Sample description:	powder - white
Sample type:	RM-reference material
Comments ¹ :	Chiron Lot#11715; Chiron
Date of entry:	8/22/2015

Substance identified-structure ² (base form)	
Systematic name:	7-bromo-5-(2-chlorophenyl)-4-hydroxy-3,4-dihydro-2H-1,4-benzodiazepin-2-one Phenazepam-4-oksido (Chiron name), Phenazepam-4-hydroxy
Other names:	
Formula (per base form)	C ₁₅ H ₁₀ BrClN ₂ O ₂
M _w (g/mol)	365.61
Salt form:	base
StdInChIKey	YNKDPCTZJPAZJC-UHFFFAOYSA-N
Compound Class	Benzodiazepines
Other active cpd. detected	phenazepam
Add.info (purity..)	99.5 % (?)

¹ This report has been produced with the financial support of the Prevention of and fight against crime Programme of the European Union (grant agreement number JUST/2013/ISEC/DRUGS/AG/6413). The contents of this report are the sole responsibility of the National Forensic Laboratory and can in no way be taken to reflect the views of the European Commission.

² Created by OPSIN free tool: <http://opsin.ch.cam.ac.uk/> DOI: 10.1021/ci100384d



Report updates

date	comments (explanation)

Supporting information

Analytical technique:	applied	remarks
GC-MS (EI ionization)	+	NFL GC-RT (min): 13.12 BP(1): 286; BP(2): 288,BP(3) :330,
FTIR-ATR	+	direct measurement
GC-IR (condensed phase)	+	

GC-MS (Agilent):

GC-method is RT locked to tetracosane (RT=9.53 min).

Injection volume 1 ml and split mode (1:50) .

Injector temperature: 280 °C.

Chromatographic separation

Column: HP1-MS (100% dimethylpolysiloxane), length 30 m, internal diameter 0.25 mm, film thickness 0.25 mm.

Carrier gas He: flow-rate 1.2 ml/min. GC oven program: 170 °C for 1 min, followed by heating up to 293 °C at a rate of 18 °C/min, hold for 6.1 min, then heating at 50 °C/min up to 325 °C and finally 2.8 min isothermal.

MSD source EI = 70 eV. GC-MS transfer line T= 235°C, source and quadrupole temperatures 280°C and 180°C, respectively. Scan range m/z scan range: from 50 (40) to 550 amu.

FTIR-ATR (Perkin Elmer): scan range 4000-400 cm⁻¹; resolution 4cm⁻¹

GC- (MS)-IR condensed phase (GC-MS (Agilent) & IR (Spectra analyses-Danny) IR scan range 4000 to 700, resolution 4cm⁻¹

GC-method:

Injection volume 1 ml and split mode (1:5) .

Injector temperature: 280 °C.

Chromatographic separation

Column: HP1-MS (100% dimethylpolysiloxane), length 30 m, internal diameter 0.25 mm, film thickness 0.25 mm.

Carrier gas He: flow-rate 1.2 ml/min. GC oven program: 170 °C for 1 min, followed by heating up to 293 °C at a rate of 18 °C/min, hold for 6.1 min, then heating at 50 °C/min up to 325 °C and finally 2.8 min isothermal.


Split MS : IR : (1:9)

MSD source EI = 70 eV. GC-MS transfer line T= 235°C, source and quadrupole temperatures 280°C and 180°C, respectively. Scan range m/z scan range: from 50 (40) to 550 amu.


IR (condensed phase): IR scan range 4000 to 700, resolution 4cm⁻¹

FIGURES OF SPECTRA

GC- MS (EI)
Abundance



FTIR-ATR



IR (condensed phase)

