

Co-funded by the Prevention of and Fight against Crime Programme of the European Union

# RESPONSE PROJECT DATABASE - NPS AND RELATED COMPOUNDS

# DESCRIPTION and GUIDELINES FOR USE

Dr. Sonja Klemenc

RS Ministry of the Interior, General Police Directorate

National Forensic Laboratory



NACIONALNI FORENZIČNI LABORATORIJ NATIONAL FORENSIC LABORATORY



Version 1.0 Ljubljana, May, 2017



Co-funded by the Prevention of and Fight against Crime Programme of the European Union

This work has been produced with the financial support of the Prevention of and fight against crime Programme of the European Union (grant agreement number JUST/2013/ISEC/DRUGS/AG/6413). The contents of this publication are the sole responsibility of the author and can in no way be taken to reflect the views of the European Commission.

# Table of contents

| 1 | Foreword   | 4  |
|---|--|----|
| 2 | Disclaimer   | 5  |
| 3 | Database description and guidelines for use                                | 6  |
|   | 3.1 Examples of use  | 10 |
| 4 | Substance descriptors  | 17 |
| 4 | 4.1 RESPONSE project structure descriptors and cross validation strategies | 19 |
| 5 | Conclusions  | 23 |
| 6 | References   | 24 |

#### 1 Foreword

The main goal of the RESPONSE<sup>1</sup> NPS part of the project was to provide and share numerous MS and FTIR spectral data on newly appearing or recently reported NPS.

Several sources of materials have been characterized and results implemented into database. Undoubtedly certified reference materials (CRM) are the most reliable source. The drawback is that CRMs are very expensive and beside this they are not available in real time, i.e. when laboratory detects a new compound. The producers are usually several steps - months behind the "gray" internet market. Thus other sources of NPS have been explored: NPS materials purchased over internet (test purchases), seized samples (police, customs), samples collected anonymously (in Slovenia) or samples collected from other laboratories.

The selection of analytical methods for chemical characterizations of samples categories listed above is dependent on the sample type and prior information available about the sample/ compound. In general the following instrumental methods have been applied: GC-MS (EI ionization 70 eV), HPLC-TOF, FTIR-ATR, GC-(MS)-IR (condensed phase), NMR (1D and 2D experiments) and later on ion chromatography (anions) and in some specific cases melting point measurements were applied as well. Interpretations of experimentally obtained data and correct identification of particular compound are checked by at minimum two experts independently. Some interpretations are supported by chem.-informatics tools. In addition, from the beginning of 2016, the structure conformation/ elucidation of the test purchased samples by NMR partner is "blind", i.e. it is done from the base. Data provided to NMR analysts are: exact monoisotopic mass, empirical formula, information of salt or base form of compound and solubility data.

Already in the first year the Slovenian partners (National Forensic Laboratory (NFL) and Faculty of chemistry and Chemical Technology (FKKT), University Ljubljana) identified many novel compounds and generated numerous MS, IR and NMR spectra. Therefore, to share gathered analytical information promptly, which means almost in real time, the public open »DRUGS MONOGRAPHS – NPS and related compounds"<sup>2</sup> was developed and lunched at the project web site at the end of 2015. Database is searchable along several criteria and at the same time it supports identification of »unknown« substance if its mass spectrum (electron ionization at 70eV) is known.

For great majority of compounds "*Analytical reports*" have been issued and implemented into database. The first part of the Analytical report describes the sample (reference material, seized, collected or test purchased sample). For reference materials specific information which guarantees traceability back to the vendor and a lot number is given. In the next session of the report analytical parameters followed by supporting experimental data obtained during chemical characterization of particular sample are given. The types and quantity of supporting experimental data depend on sample type and previous information already available from other resources (databases and/or scientific publications). In 2016 database was enriched with several reports from the contributors outside the RESPONSE project. We kindly acknowledge this and at the same time encourage the new contributors to participate!

Dr. Sonja Klemenc RESPONSE project manager

#### 2 Disclaimer

Although "RESPONSE" analytical results are reviewed independently by at minimum two experts of National Forensic Laboratory (NFL) chemistry department and NMR interpretations by two experts from Faculty of Chemistry and Chemical Technology, University of Ljubljana, it is still sole responsibility of database user to critically re-evaluate our findings. In the case of the observed inconsistencies, i.e. possible identification errors please report them to sonja.klemenc@policija.si or nfl@policija.si .

RESPONSE database shall be regarded as supporting tool only and not as an officially peer-reviewed database.

#### 3 Database description and guidelines for use

Database window is shown on Figure 1. Database is structured as a table. Each row is a record which describes one compound. Heading row currently has 18 fields – compound/sample descriptors - see details in the Table 1). Each field has a *sort button* (example is labeled by red circle no 1 on Figure 1) which enables sorting of entries (alphabetically or ascending/descending order).

On the top left side of the window there is a field *"show"* (red circle no. 2 on Figure 1). It serves for the adjustments of database window, i.e. how many records one wish to see per page. Options are: 5, 10, 25, 50 and All. To see the scroll bar at the bottom of the window we suggest the choice of option 5 or 10.

At the bottom (left corner) there is information of what you see on the screen (example on Figure 1: *Showing 1 to 5 of 422 entries*), a date of the last database update is given bellow and counter of database "visits" is at the bottom.

Buttons on the bottom right side of the window enables movements between the pages. Scroll bar enables left – right movements along the database window. Vertical scroll bar will appear on as needed basis.

| ow 5 🔹 entrie   | ~                              | 1  |  |  |                          |                                     |             |             |             |                     |                             |                                  |                          | $\left( \right)$                  | Search:          |                   | $\sum$             |
|-----------------|--------------------------------|--|--|--|--------------------------|-------------------------------------|-------------|-------------|-------------|---------------------|-----------------------------|----------------------------------|--------------------------|-----------------------------------|------------------|-------------------|--------------------|
| Substance class | substance (NPS)<br>common name | structure<br>(created by<br>OPSIN free tool) | NPS1 systematic name                                 | other names 🍦                            | Formula per<br>base form | Mw (g/mol)<br>per base<br>form NPS1 | MS<br>(BP1) | MS<br>(BP2) | MS<br>(BP3) | GC-MS-RT<br>NFL/min | MS<br>spectrum<br>(picture) | Stdin<br>ChiKey                  | Type of detection        | comments                          | date of<br>entry | report<br>updated | Chem-ID<br>Report- |
| Aminoindanes    | <u>N-Methyl-2Al</u>            | HCC)   | N-methyl-2,3-dihydro-<br>1H-inden-2-amine            | NM-2AI, N-<br>methyl-2-<br>aminoindane   | C10H13N                  | 147.22                              | 147         | 42          | 106         | 2.73                | <u>show</u>                 | SXV/ZQUCTTOBHJT-<br>UHFFFAOYSA-N | test purchase            |                                   | 2015-<br>10-25   | 2015-12-<br>03    | 1349-15            |
| Aminoindanes    | MDAi                           | <u> </u>                                     | 2H,5H,6H,7H-<br>indeno[5,6-d]<br>[1,3]dioxol-6-amine | 5,6-<br>methylenedioxy-<br>2-aminoindane | C10H11NO2                | 177.20                              | 160         | 149         | 177         | 4.31                | <u>show</u>                 | FODRMHHCWZAXJM-<br>UHFFFAOYSA-N  | test purchase            |                                   | 2015-<br>11-25   | 2015-12-<br>03    | 1347-15            |
| Aminoindanes    | <u>2-Al</u>                    | ())  | 2,3-dihydro-1H-inden-<br>2-ylamine                   | 2-Aminoindane                            | C9H11N                   | 133.19                              | 133         | 116         | 91          | 2.55                | <u>show</u>                 | LMHHFZAXSANGGM-<br>UHFFFAOYSA-N  | RM-reference<br>material | I-SEE<br>project RM<br>and report | 2016-<br>01-20   |                   | 1469-16            |
| Aminoindanes    | <u>5-IAI</u>                   | <u>`</u> CD-                                 | 5-lodo-2,3-dihydro-1H-<br>inden-2-amine              | 5-lodo-2-<br>aminoindane                 | C9H10IN                  | 259.09                              | 259         | 115         | 117         | 4.86                | <u>show</u>                 | BIHPYCDDPGNWQO-<br>UHFFFAOYSA-N  | RM-reference<br>material | I-SEE<br>project RM<br>and report | 2016-<br>01-20   |                   | 1475-16            |
| Aminoindanes    | <u>2-AI</u>                    | ()))   | 2,3-dihydro-1H-inden-<br>2-amine                     | 2-Aminoindane                            | C9H11N                   | 133.19                              | 133         | 116         | 91          | 2.50                | <u>show</u>                 | UMHHFZAXSANGGM-<br>UHFFFAOYSA-N  | RM-reference<br>material |                                   | 2016-<br>07-18   |                   | 1598-16            |
|                 |                                |  |  |  |                          |                                     |             |             |             |                     |                             |                                  |                          |                                   |                  |                   |                    |

Figure 1: Database window at 28/10/2016

On the top right side of the window (red circle no. 3 on Figure 1) there is a field *"search*". You can search for any reasonable parameter. The search field functions as a filter.

|   | Field                                   | description  | Comments/actions-use  | notes  |
|---|---|--|---|--|
| 1 | Substance class                         | Substances are classified into several classes along their physicochemical characteristics /aminoindanes, arylalkylamines, arylcyclohexylamines, benzodiazepines, cannabinoids, cathinones, indolalkylamines (tryptamines), opioids, others, plants & extracts.<br>This classification was adopted from EMCDDA (in 2015). <sup>3</sup> | Click the <b>sort button</b> – records<br>will be sorted from A-Z or in<br>opposite direction;<br>By typing the name of particular<br>class into the <b>search field</b> you<br>can select all records from<br>particular class and then use<br>the <b>sort buttons</b> in other fields<br>for deeper search. |  |
| 2 | Substance<br>Common name                | Be aware that this is not a unique<br>name (for the same substance we can<br>find more than one common name in<br>different data repositories). For test<br>purchased compounds we (NFL often<br>follow the name given by web wendor.  | Click the <b>sort button</b> – records<br>will be sorted from A-Z or in<br>opposite direction;<br>Click on the link (compound<br>name) to see the <b>"Analytical</b><br><b>report"</b> (pdf format)   |  |
| 3 | Structure                               | Picture of structure (in .png format)  | Click or double click to enlarge<br>– structure opens in new<br>window.   | Structure – name – structure was<br>crosschecked by free tools available<br>over the internet (procedures are<br>explained later on).  |
| 4 | NPS1<br>systematic name                 | Systematic name – IUPAC based<br>name which describes the compound<br>structure.   | Click the <b>sort button</b> to sort<br>records ascending/ descending<br>order.<br>Apply the search field to search<br>for the whole name or only a<br>part of the name (for example<br>indole or indazole or<br>carboxylate)   | In the RESPONSE project data<br>systematic name is always given<br>per base form of a compound.<br>Structure – name –structure is<br>crosschecked by free tools available<br>over the internet (procedures are<br>explained later on). |
| 5 | Other names                             | Other names (or some of them) for<br>compound which are available over<br>the internet at the time of the<br>compound entry (IUPAC based or<br>common names).  | Click the <b>sort button</b> – to sort<br>records ascending/ descending<br>order.<br>Apply the search field to search<br>for the whole or only a part of<br>the name (for example AKB.or<br>PICA or PINACA, etc).   |  |
| 6 | Formula per<br>base form                | Empirical formula of the compound per base form.   | Click the <b>sort button</b> to sort the records (ascending/descending) along the number of C atoms.  |  |
| 7 | M <sub>w</sub> per base<br>form (g/mol) | Molecular weight per base form.  | Click the <b>sort button</b> to sort the records ascending/descending).   |  |

# Table 1: Database fields description and guidelines for use

|    | Field                    | description   | Comments/actions-use  | notes   |
|----|--------------------------|---|---|---|
| 8  | GC-MS-RT<br>NFL (min)    | Retention time (RT) as measured in the National forensic laboratory (NFL).  | Click the <b>sort button</b> to sort the records.   | In NFL we use the "retention time<br>locked" GC-MS method. Therefore,<br>absolute RT shows long term<br>reproducibility at our analytical<br>conditions. For other user this field<br>is not very useful. |
| 9  | MS (BP1)                 | The most intensive peak (m/z) in the mass spectrum.   | Click the <b>sort button</b> to sort the records (ascending/ descending order).   | For NFL data the peak areas calculated automatically (by Agilent Software) are taken into account.  |
| 10 | MS (BP2)                 | Second most intensive peak (m/z) in the mass spectrum.  | Click the <b>sort button</b> to sort the records (ascending/ descending order)  | Be aware that isotopic peaks are taken into account as well.  |
| 11 | MS (BP3)                 | Third most intensive peak (m/z) in the mass spectrum.   | Click the <b>sort button</b> to sort the records (ascending/ descending order)  |   |
| 12 | MS spectrum<br>(picture) | Mass spectrum (E, +70eV) picture (in .png format).  | Click to see the picture which opens in new window  | Majority of MS spectra were<br>collected in NFL. Exception are the<br>spectra provided from contributors<br>of analytical reports (see type of<br>detection field).                                       |
| 13 | StdInChlKey              | Standard InChI Key = Standardized<br>International Chemical Identifier is a<br>textual identifier for small organic<br>molecules. The condensed, 27<br>characters Standard InChIKey is a<br>hashed version of the full standard<br>InChI string, designed to allow for<br>easy web searches of chemical<br>compounds. <sup>4</sup>  | Click <b>the link</b> and search the<br>internet for information on<br>particular compound.   | Code is crosschecked by free tools<br>available over the internet. See<br>more about coding in the heading<br>4.1).   |
| 14 | Type of<br>detection     | The source of compound (sample) or<br>analytical information<br>sources:<br><u>Analytical report</u> – analytical data were<br>provided from contributors outside<br>NFL. For these records NFL does not<br>have the sample.<br><u>Collected</u> – samples collected<br>anonymously in Slovenia or provided<br>from other laboratories<br><u>In house RM</u> – in house reference<br>materials obtained in the frame of<br>proficiency tests or assessment of the<br>lab<br><u>RM-reference materials</u> – purchased<br>from official vendors<br><u>Seized</u> – samples seized by the<br>Slovenian Police.<br><u>Test purchase</u> – purchased samples<br>which are advertised and sold over the<br>internet as legal highs or research<br>chemicals (not reference materials!) | Click the sort button to sort the<br>records alphabetically;<br>By typing the name of particular<br>"type of detection" into the<br>search field you will see all<br>records of the same type of<br>detection. Afterwards you can<br>use the sort buttons in other<br>fields for deeper search. | Analytical strategy is dependent on<br>the type of detection and is<br>described in a separate document.<br>(see below)   |
| 15 | Comments                 | Comments about the analytical findings and/or contributor of the sample or analytical report  | Click the sort button to sort ascending/descending.   |   |
| 16 | Date of entry            | Date when the record was entered into<br>a database   | Click the sort button to sort ascending/descending  |   |

|    | Field                   | description   | Comments/actions-use                                  | notes   |
|----|-------------------------|---|---|---|
| 17 | Report updated          | Date of the "Analytical report" update  | Click the sort button to sort ascending/descending    | Some "cosmetic" improvements of the reports will not be labeled as updates.   |
| 18 | Chem-ID or<br>Report-ID | Chem-ID is the unique label of<br>particular sample in the NFL chemicals<br>database.<br>REPORT-ID is the unique label of the<br>analytical report provided by<br>contributors outside NFL. | Click the sort button to sort<br>ascending/descending | This code guarantees the traceability to the sample and/or original source of chromatographic, spectroscopic and other data provider for the particular sample. |

#### 3.1 Examples of use

#### Example 1: The SEARCH field

One is interested only on reference materials records, sorted by molecular weight.

Procedure:

- in the "search" field type "reference"
- apply "sort" button in M<sub>w</sub> (g/mol) field

×

Database respond as obtained at 31/10/2016: 222 entries (filtered from 422 total entries), sorted in ascending values of  $M_w$  (Figure 2).

To remove the current filter click on blue cross in the search field or apply the F5 key on the keyboard.

Search: reference

| w 5 • entries      |                                |  |  |   |                          | 0                                   |             |             |             |                     |                             |                                  |                          | Search: refer                  | ence             |
|--------------------|--------------------------------|--|--|---|--------------------------|-------------------------------------|-------------|-------------|-------------|---------------------|-----------------------------|----------------------------------|--------------------------|--------------------------------|------------------|
| Substance<br>class | substance (NPS)<br>common name | structure<br>(created by<br>OPSIN free tool) | NPS1 systematic name                       | other names   | Formula per<br>base form | Mw (g/mol)<br>per base<br>form NPS1 | MS<br>(8P1) | MS<br>(BP2) | MS<br>(BP3) | GC-MS-RT<br>NFL/min | MS<br>spectrum<br>(picture) | Stdin<br>ChiKey                  | Type of detection        | comments                       | date of<br>entry |
| thers              | Indane                         | (1)  | 2,3-dillydro-1H-indene                     | 1,2-hydrindene,<br>hydrindonaphthe<br>ne,<br>hydrindene                               | C9H10                    | 118.18                              | 117         | 118         | 115         | 2.09                | show                        | PONFLJBRINBOBRO:<br>UHFFFAOYSA-N | RM-reference<br>material |                                | 2016-08-03       |
| minoindanes        | <u>2-AI</u>                    | 00-  | 2,3-dihydro-1H-inden-<br>2-ylamine         | 2-Aminoindane   | C9H11N                   | 133,19                              | 133         | 116         | 91          | 2.55                | stow                        | LWHHEZAXSANGGM-<br>UHEFEAOYSA-N  | RM-reference<br>material | I-SEE project RM<br>and report | 2016-01-20       |
| minoindanes        | <u>2-AI</u>                    | 00-  | 2,3-dihydro-1H-inden-<br>2-amine           | 2-AminoIndane   | C9HTIN                   | 133.19                              | 133         | 116         | 91          | 2.50                | show                        | LMHHEZAXSANGGM-<br>UHEEFAQYSA-N  | RM-reference<br>material |                                | 2016-07-18       |
| thers              | 1,3-Indanedione                | B  | 2,3-dihydro-1H-indene-<br>1,3-dione        | 1,3-Dicxolndan  | C9H6O2                   | 146.15                              | 146         | 104         | 76          | 2.90                | stow                        | UHKAJLSKXRADET-<br>UHEEFAOYSA-N  | RM-reference<br>material |                                | 2016-08-03       |
| minoindanes        | N-methyl-2Al                   | 400  | N-methyl-2,3-dihydro-<br>1H-inden-2-antine | 2,3-dihydro-N-<br>methyl-1H-<br>inden-2-amine<br>NW-2AI<br>N-methyl-2-<br>Aminoindane | C10H13N                  | 147.22                              | 147         | 42          | 106         | 2.75                | show                        | SXVZQUCTTORHJT-<br>UHFFFAOYSA-N  | RM-reference<br>material |                                | 2016-07-18       |

Figure 2: Filtered records on reference materials and sorted by M<sub>w</sub>.

**Example 2:** Identify your "unknown" by searching RESPONSE project database MS spectral data (BP1, BP2, BP3) and picture of MS spectrum. This example comes from the real request.

Let's assume that a expert has the mass spectrum as shown on Figure 3 in hand. Mw is not known, the compound class is not known.

#### First step:

Read peak intensities from your spectrum: BP1=144, BP2=144, BP3=233




Figure 3: Mass spectrum of the "unknown" compound (real request)

#### Second step:

Search the RESPONSE project database:

- 3.1.1 Clear all previous selections from the "search field"
  - a. In the "search field" type: 232 144 233 (put spaces between numbers)




Figure 4: Search field criteria

3.1.2 Result (as obtained at 20/11/2016) is shown in Figure 5: *Six compounds* had been filtered out of 422 entries

| Substance<br>class | substance (NPS)<br>common name        | structure<br>(created by<br>OPSIN free tool) | NPS1 systematic name   | other names  | Formula per<br>base form | Mw (g/mol)<br>per base<br>form NPS1 | MS<br>(BP1) | MS<br>(BP2) | MS<br>(BP3) | GC-MS-RT<br>NFL/min | MS<br>spectrum<br>(picture) | Stdin<br>ChiKey                                |
|--------------------|---------------------------------------|--|--|--|--------------------------|-------------------------------------|-------------|-------------|-------------|---------------------|-----------------------------|--|
| Cannabinoids       | 5F-MDMB-PICA                          | 18th   | methyl 2-{[1-(5-<br>fluoropentyl)-1H-indol-<br>3-yl]formamido]-3,3-<br>dimethylbutanoate | 5F-MDMB-2201,<br>MDMB-2201   | C21H29FN2O3              | 376.46                              | 232         | 144         | 233         |                     | <u>show</u>                 | <u>CHSUEEBESACQDV-</u><br>UHEFFAOYSA-N         |
| Cannabinoids       | <u>5F_PB-22</u>                       | Sec.   | quinolin-8-yl 1-(5-<br>fluoropentyl)-1H-<br>indole-3-carboxylate                         | 5-fluoro QUPIC,<br>5-fluoro-PB-22  | C23H21FN202              | 376.43                              | 232         | 144         | 233         | 17.11               | <u>show</u>                 | MBOCMBEDYVSGLJ-<br>UHFFFAOYSA-N                |
| Cannabinoids       | 5-fluoro-JVH-018-<br>adamantyl analog | \$*\$ر_                                      | 3-(adamantane-1-<br>carbonyl)-1-{5-<br>fluoropentyl)-1H-<br>indole                       | [1-(5-<br>fluoropentyl)-1H<br>-indol-3-<br>yl]tricyclo[3.3.1.<br>13,7]dec-1-yl-<br>methanone;<br>5F-AB-001;<br>AM2201<br>adamantyl<br>analog | C24H30FNO                | 367.51                              | 232         | 144         | 233         | 15.37               | <u>show</u>                 | <u>WBLLMFXPGHZFFP-</u><br>UHFFFAOYSA- <u>N</u> |
| Cannabinoids       | <u>LTI-701</u>                        | 2840   | 1-(5-fluoropentyl)-N-<br>phenyl-1H-indole-3-<br>carboxamide                              |  | C20H21FN2O               | 324.40                              | 232         | 144         | 233         | 14.05               | <u>show</u>                 | FZYJALPWLNXEIP-<br>UHFFFAOYSA-N                |
| Cannabinoids       | <u>LTI-701</u>                        | 1.840  | 1-(5-fluoropentyl)-N-<br>phenyl-1H-indole-3-<br>carboxamide                              |  | C20H21FN2O               | 324.40                              | 232         | 144         | 233         | 14.01               | <u>show</u>                 | <u>FZYJALPWLNXEIP-</u><br><u>UHFFFAOYSA-N</u>  |

Showing 1 to 6 of 6 entries (filtered from 422 total entries) Table updated 28.10.2016

19748 hits

Figure 5: Filtered entries for the criteria given in search field (remark: only five records are shown on the figure above).

- b. To inspect the mass spectra of filtered compounds and compare them with your data click the "show" button in the column "MS spectrum picture" (picture opens in new window).
- c. RESULT: For this particular example the mass spectrum of 5F-MDMB-PICA (Figure 6) was in good agreement with the spectrum in question.




Figure 6: 5F-MDMB-PICA spectrum in the RESPONSE project database

#### Third step:

See »Analytical report« (click the compound common name link) to see if some other analytical data are presented (as for example FTIR and other analytical data).

To search the WEB for information on this substance – click the StdInChI Key link.

| Substance<br>class | substance (NPS)<br>common name         | structure<br>(created by<br>OPSIN free tool) | NPS1 systematic name   | other names 🍦  | Formula per<br>base form | Mw (g/mol)<br>per base<br>form NPS1 | MS<br>(BP1) | MS<br>(BP2) | MS<br>(BP3) | GC-MS-RT<br>NFL/min | MS<br>spectrum<br>(picture) | Stdln<br>ChlKey                               |
|--------------------|--|--|--|--|--------------------------|-------------------------------------|-------------|-------------|-------------|---------------------|-----------------------------|---|
| Cannabinoids       | 5F-MDMB-PICA                           | 18th   | methyl 2-{[1-(5-<br>fluoropentyl)-1H-indol-<br>3-yl]formamido]-3,3-<br>dimethylbutanoate | 5F-MDMB-2201,<br>MDMB-2201   | C21H29FN2O3              | 376.46                              | 232         | 144         | 233         | )                   | show                        | CHSUEEBESACODV-<br>UHFFFAOYSA-N               |
| Cannabinoids       | <u>5F_PB-22</u>                        | SS .   | quinolin-8-yl 1-(5-<br>fluoropentyl)-1H-<br>indole-3-carboxylate                         | 5-fluoro QUPIC,<br>5-fluoro-PB-22  | C23H21FN202              | 376.43                              | 232         | 144         | 233         | <b>17.</b> 11       | <u>show</u>                 | MBOCMBFDYVSGLJ-<br>UHFFFAOYSA-N               |
| Cannabinoids       | 5-fluoro-JV/H-018-<br>adamantyl analog | <u>,</u> \$*\$                               | 3-(adamantane-1-<br>carbonyl)-1-{5-<br>fluoropentyl}-1H-<br>indole                       | [1-(5-<br>fluoropentyl)-1H<br>-indol-3-<br>yl]tricyclo[3.3.1.<br>13,7]dec-1-yl-<br>methanone;<br>5F-AB-001;<br>AM2201<br>adamantyl<br>analog | C24H30FNO                | 367.51                              | 232         | 144         | 233         | 15.37               | <u>show</u>                 | WBLLMFXPGHZFFP-<br>UHFFFAQYSA-N               |
| Cannabinoids       | <u>LTI-701</u>                         | 2800   | 1-(5-fluoropentyl)-N-<br>phenyl-1H-indole-3-<br>carboxamide                              |  | C20H21FN2O               | 324.40                              | 232         | 144         | 233         | 14.05               | <u>show</u>                 | FZYJALPWLNXEIP-<br>UHFFFAOYSA-N               |
| Cannabinoids       | <u>LTI-701</u>                         | 'Sto   | 1-(5-fluoropentyl)-N-<br>phenyl-1H-indole-3-<br>carboxamide                              |  | C20H21FN2O               | 324.40                              | 232         | 144         | 233         | 14.01               | <u>show</u>                 | <u>FZYJALPWLNXEIP-</u><br><u>UHFFFAOYSA-N</u> |

Remark: If the molecular weight and/or empirical formula are known the number of possible candidates can be reduced. For molecular weight search we recommend to use the whole number (neglect decimal digits).

**Example 3:** You have formal name of a compound as for example: 1-[(tetrahydro-2H-pyran-4-yl)methyl]-N-tricyclo[3.3.1.13,7]dec-1-yl-1H-indazole-3-carboxamide

Does the compound exist in the RESPONSE database?

<u>First step</u>: Apply the "search field" and search for: 1-[(tetrahydro-2H-pyran-4-yl)methyl]-N-tricyclo[3.3.1.13,7]dec-1-yl-1H-indazole-3-carboxamide

There is no match in the database. Is this correct? Go to the next step and try another way around.

<u>Second step</u>: Go to the internet, apply OPSIN application (<u>http://opsin.ch.cam.ac.uk/</u>) and generate the StdInChI Key from the IUPAC name.




Figure 7: StdInChIKey from formal name

Third step: Copy the generated code (*WFDRVDUEMGJILX-UHFFFAOYSA-N*) into the RESPONSE database "Search field".

Result: ADAMANTYL-THPINACA (see Figure 8).

| Show All 🔻 entries   |                                 |  |  |             |                          |                                       |             |             |             |                     |                             |                                 |
|--|---------------------------------|--|--|-------------|--------------------------|---------------------------------------|-------------|-------------|-------------|---------------------|-----------------------------|---------------------------------|
| Substance<br>class   | substance (NPS)<br>common name  | structure<br>(created by ∲<br>OPSIN free tool) | NPS1 systematic name   | other names | Formula per<br>base form | Mw (g/mol)<br>per base ≑<br>form NPS1 | MS<br>(BP1) | MS<br>(BP2) | MS<br>(BP3) | GC-MS-RT<br>NFL/min | MS<br>spectrum<br>(picture) | Stdin<br>ChiKey                 |
| Cannabinoids   | ADAMANTYL-<br>THPINACA          | (D) rg?  | N-(1-adamantyl)-1-<br>(tetrahydropyran-4-<br>ylmethyl)indazole-3-<br>carboxamide |             | C24H31N3O2               | 393.52                                | 135         | 243         | 145         | 17.52               | <u>show</u>                 | WFDRVDUEMGJILX-<br>UHFFFAOYSA-N |
| Showing 1 to 1 of 1 entr<br>Table updated 28.10.20<br>19748 hits | ies (filtered from 422 to<br>16 | tal entries)                                   |  |             |                          |                                       |             |             |             |                     |                             |                                 |

Figure 8: Search (at 20/11/2016) along StdInChIKey returns compound ADAMANTYL THPINACA.

#### 4 Substance descriptors

Each distinct compound has a unique molecular structure which can be designated by a structural formula. Therefore, the expectation is that each structure will be "accompanied" by characteristic and unique name and/or other unique descriptor recognizable by humans and by computerized systems and different types of electronic databases, either local or public opened. The unique descriptor – substance identifier enables consistency of substance coding within and between databases.

Common names and synonyms are not unique and mostly not related to the compound structure. Therefore, they are "poor descriptors" of the structure but widely used in conversations. It is known from literature<sup>5</sup> that some authority-assigned chemical identifiers, which are produced by assignment made by database keepers like Chemical Abstracts Service (CAS) Registry numbers by the American Chemical Society<sup>6</sup>, EC numbers from the European Community Inventory<sup>7</sup>, CID and SID numbers assigned by PubChem<sup>8</sup>, ChemSpider<sup>9</sup>, etc are not the best solution. Namely, different databases use different codes for the same structure and for many novel NPS such codes do not exist at all. Well, some converters<sup>10</sup> can be helpful for finding links of between databases encoding systems.

The alternative to authority-assigned codes and common names are structure-based chemical identifiers. These are derived from structural formula. Anyhow, systematic names (IUPAC based) generated for the same structure often differ. SMILES (simplified molecular-input line-entry system)<sup>11</sup> generated for the same structure by different encoding tools shows great diversity. A review of these and other structure notations is given in work of Warr<sup>12</sup>. It seems that currently the best way to annotate the structure uniquely is it's encoding by the International Chemical Identifier code, i. e. InChI string or by its hashed vesion InChI Key. In-depth discussion, some advantages and/or disadvantages one can find in reference<sup>5</sup> and references listed there.




Figure 9: Standard InChI string and corresponding standard InChI Key – example for methcathinone (2-(methylamino)-1-phenyl-propan-1-one).

# 4.1 **RESPONSE** project structure descriptors and cross validation strategies

During the RESPONSE project database development and further evolution we explored several types of structure descriptors<sup>13</sup>. Namely, the idea was that each structure should be:

Described by several generally accepted descriptors in the NPS field investigations

- structure (picture as primary descriptor)
- structure related descriptor IUPAC name (English only)
- structure encoding that machines and humans can understand and/or compare i.e.
  systematic chemical identifier should be given and it should be "unique" descriptor ,
  if possible
- o other descriptors (common names and/or synonyms)...

The "unique" descriptor should:

- o enhance interoperability options with other public open databases
- be "googlable" (it can serve as self-updating info tool i. e. new info about particular compound can be reached in real time)

Validation of key identifiers (structure-name-structure) should be performed by Chem-informatic tools.

4.1.1 Unique descriptor


At the beginning RESPONSE used canonical SMILES, which are not unique. Because of this SMILES were replaced by "Standard InChI Key", which in our opinion is good approximation to the "unique" descriptor.

The main advantage of InChI Key over the InChI string is its uniform length (27 characters – see Figure 9) for any structure, while the string code is longer, of flexible length and contains some symbols not suited the best for search engines and web browsers.

The same InChI label always means the same substance, and the same substance always receives the same label, under the same labeling conditions – <u>for the same level of structure details</u>!

Figure 10 shows the same compound with different levels of structure details (a) and (b) which gives different codes. Difference is the second (8 letters block), which encodes stereochemistry and

isotopes. The first block of letters, which defines skeleton (connectivity) is identical in both examples. Therefore, one can apply the first block of letters (ROWZIXRLVUOMCJ) in search engines (google for example) to see information on compound regardless of its stereochemistry. Try it!



Structure b): <u>ROWZIXRLVUOMCJ-FQEVSTJZSA-N</u>

The original InChI trust software (free software and free public open source code) was applied for Standard InChI key generation from the structure (sdf or mol format) and OPSINE free tool for code generation from the IUPAC name.

In the RESPONSE project database Standard InChI Key is given:

- for the **base form** of compound
- for un-charged state of compound
- for test purchased, collected and seized samples the stereochemistry is not taken into account (i.e. it is undefined or unknown), while for reference materials data from certificates are followed (strictly from the end of 2016).

One should be aware that codes for base and salt forms, stereoisomers, tautomers and charged species differ!

Figure 10: The same compound with different levels of structure details (a) and (b) gives different codes. Difference is the second (8 letters blok), which encodes stereochemistry and isotopes. The first block of letters, which defines skeleton connectivity is identical in both examples.

#### 4.1.2 Validation of key identifiers

In the RESPONSE project we explored several free tools for the conversion of structure (in the case of seized, collected and test purchased samples) into IUPAC based name or opposite (when reference materials were in question) and for generation of standard InChIKey code and for cross-validation validation of these data as well. The cross validation process is schematically depicted in figure below.




Figure 11: Schematic presentation of RESPONSE project cross-validation tools explored. (blue shadowed parts of the scheme were omitted from the finally implemented validation strategy)

<u>Marvin Sketch</u><sup>14</sup> (ChemAxon) was our basic tool for structural formula generation followed by saving the structure as a electronic "molfile" (readable by other applications), generation of systematic IUPAC and/or traditional name, Mw, exact mass and some other calculations. In addition for cross validations we explored and combined the following free tools:

<u>INCHI-TRUST</u><sup>15</sup> standard and non-standard InChI/InChIKey generator (free software) which is available here: <u>http://www.inchi-trust.org/downloads/</u>. This software generates codes from **structural formula** (input is .mol file), which, in our opinion, is the correct way of code generation. So far, we applied the tool on more than 400 entries for cross validation of standard InChIKey obtained from the OPSIN application. We did not find any inconsistencies. Anyhow, the advantage of INCHI-TRUST software is that user can define InChI processing options<sup>16</sup>.

<u>OPSIN</u><sup>17</sup> (Open Praser for Systematic IUPAC Nomenclature, Cambridge University) was applied to generate a structure and a standard InChI/InChIKey from a systematic name. Structure obtained from OPSIN was compared with the one from Marvin Sketch. Standard InChIKey was validated by the one generated from a structure by INCHI-TRUST generator. Picture of structure in ".png" format and generated StdInChIKey were implemented into RESPONSE database fields.

<u>CACTUS</u> (NCI/CADD Group)<sup>18</sup> Chemical identifier resolver, which allows the conversion of a given structure identifier into another representation or structure identifier is available here:

<u>https://cactus.nci.nih.gov/chemical/structure</u>. This application was applied for cross validation of Standard InChI Key generated by OPSIN. In some rear cases we observed inconsistencies between Std InChiKey gererated

<u>CHEMICALIZE<sup>19</sup></u> (ChemAxon) is a powerful online platform for chemical calculations, search and text processing. In the "Calculate" module it generates structure from name and/or name from structure, InChI/InChIKey and performs many other calculations: <u>https://chemicalize.com/#/calculation</u>. We applied this toll for different purposes. However, with regard to InChI key we observed a "drawback". Namely, in some occasions software generates standard InChI string and/or Key and sometimes nonstandard InChI Key. Along the best of our knowledge user has no influence on this. From the encoding result it can be deduced which type of identifier was actually generated (see Figure 12). Anyhow, for consistency within and between databases such "random" code generation does not fit. This was the main reason why we ruled out this tool (or MarvinSketch directly) for InChI Key generation.

| Names and identifie | rs A  |
|---------------------|---|
| Common names        | methcathinone   |
| IUPAC name          | 2-(methylamino)-1-phenylpropan-1-one                          |
| SMILES              | CNC(C)C(=O)C1=CC=CC=C1  |
| InChl               | InChI=1/C10H13NO/c1-8(11-2)10(12)9-6-4-3-5-7-9/h3-8,11H,1-2H3 |
| InChlKey            | InChIKey=LPLLVINFLBSFRP-UHFFFAOYNA-N                          |
| CAS                 | 28521-94-0,5650-44-2  |

| Names and identifiers | В  |
|-----------------------|--|
| Common names          | N/A  |
| IUPAC name            | 1-pentyl-N-phenyl-1H-indole-3-carboxamide  |
| SMILES                | CCCCCN1C=C(C(=O)NC2=CC=CC2)C2=CC=CC=C12  |
| InChl                 | InChI=15/C20H22N2O/c1-2-3-9-14-22-15-18(17-12-7-8-13-19(17)22)20(23)21-16-10-5-4-6-11-16/h4-8,10-13,15H,2-3,9,14H2,1H3,(H,21,23) |
| InChlKey              | InChIKey=SPUQBSMPSCYQLS-UHFFFA0YSA-N   |
| CAS                   | N/A  |

Figure 12: Chemicalize outputs for two compounds: A) for methcathinone (2-(methylamino)-1-phenyl-propan-1-one) **non standard** InChl string and Key were generated; B9) for SDB-006-N-phenyl-analog (1-pentyl-N-phenyl-1H-indole-3-carboxamide) **standard** forms of string and Key were the output.

### 5 Conclusions

The RESPONSE project database tool is technically relatively simple and far from being perfect. It was created ad hoc and at no costs. Anyhow, it's "designers" hope that it fits for the intended purpose of use.

If you have found RESPONSE database "NPS AND RELATED COMPOUNDS" useful in your work citing our database would be very much appreciated.

#### 6 References

<sup>5</sup> Stephen R Heller, Alan McNaught, Igor Pletnev, Stephen Stein and Dmitrii Tchekhovskoi, InChI, the IUPAC International Chemical Identifier, Journal of Cheminformatics 2015, 7:23, DOI: 10.1186/s13321-015-0068-4; <a href="https://jcheminf.springeropen.com/articles/10.1186/s13321-015-0068-4">https://jcheminf.springeropen.com/articles/10.1186/s13321-015-0068-4</a>

<sup>6</sup> CAS REGISTRY and CAS RN FAQs. <u>http://www.cas.org/content/chemical-substances/faqs</u>

<sup>7</sup> EC Inventory - ECHA. <u>http://echa.europa.eu/in/information-on-chemicals/ec-inventory</u>

<sup>8</sup> https://pubchem.ncbi.nlm.nih.gov/help.html

<sup>9</sup> ChemSpider | Search and share chemistry. <u>http://www.chemspider.com/</u>

<sup>10</sup> <u>http://cts.fiehnlab.ucdavis.edu/</u>

<sup>11</sup> https://en.wikipedia.org/wiki/Simplified\_molecular-input\_line-entry\_system

<sup>12</sup> Wendy Warr, Representation of chemical structures. Wiley Interdiscip Rev Comput Mol Sci. 2011;1:557–79, DOI: 10.1002/wcms.36

<sup>13</sup> S. Klemenc, Substance descriptors: structure, common name, systematic name, systematic chemical identifiers and their consistency, (presented by D.Saboti), ENFSI-DWG meeting 2017, Linköping, Sweden <a href="http://www.policija.si/eng/images/stories/GPUNFL/PDF/2017/ENFSI-">http://www.policija.si/eng/images/stories/GPUNFL/PDF/2017/ENFSI-</a>

DWG\_Meeting2017\_SubstancDescriptor.pdf

<sup>14</sup>https://www.chemaxon.com/products/marvin/marvinsketch/?gclid=CjwKEAiAgavBBRCA7ZbggrLSkUcSJACW DexAN5B9vnh01oHXQKgvgUvcMGZY065FgDNVurKrfl9VWxoC-sjw\_wcB

<sup>15</sup> http://www.inchi-trust.org/

<sup>16</sup> IUPAC International Chemical Identifier (InChI) Programs InChI version 1, software version 1.04 (September 2011) User's Guide, download able here: <u>http://www.inchi-trust.org/downloads/</u> (*INCHI-1-DOC.zip*)

<sup>17</sup> http://opsin.ch.cam.ac.uk/

<sup>18</sup> https://cactus.nci.nih.gov/

<sup>19</sup> <u>https://chemicalize.com/welcome</u>

<sup>&</sup>lt;sup>1</sup> <u>http://www.policija.si/eng/index.php/generalpolicedirectorate/1669-nfl-page-response</u>

<sup>&</sup>lt;sup>2</sup> http://www.policija.si/apps/nfl\_response\_web/seznam.php

<sup>&</sup>lt;sup>3</sup> http://www.emcdda.europa.eu/

<sup>&</sup>lt;sup>4</sup> https://en.wikipedia.org/wiki/International\_Chemical\_Identifier